

Operation Manual
for
GAS SELECT MODULE GSM (73-1647)
and
PRESSURE REGULATOR PRM (73-0195)

(Version 1.2 / TB / Printed: 07 / 2012)



1. Introduction, manufacturer's details

These Operating Instructions describe the function and use of the Gas Select Module GSM and Pressure Regulator Module PRM. They form part of the set have to be kept close to it.

All the information in these Instructions has been drawn up after careful examination but does not represent a warranty of product properties. Alterations in line with technical progress are reserved.

This device manufactured by:

**Hugo Sachs Elektronik -
HARVARD APPARATUS GmbH**
Gruenstrasse 1
79232 March- Hugstetten / Germany

Phone (Germany) : 07665-9200-0
Phone (others) : (int + 49) 7665-9200-0
Fax (Germany) : 07665-9200-90
Fax (others) : (int + 49) 7665-9200-90
eMail: sales@hugo-sachs.de

1.1 Copyright

This product and the corresponding documentation are protected by copyright. All rights reserved. This document must not be copied, photocopied, reproduced or translated, either as a whole or in parts, without prior written agreement by HUGO SACHS ELEKTRONIK - HARVARD APPARATUS GmbH, March/Hugstetten, Germany.

1.2 Trademarks

PLUGSYS® is a registered trademark of HUGO SACHS ELEKTRONIK - HARVARD APPARATUS GmbH, March-Hugstetten. Other trademarks shown in the Operating Instructions are the property of the corresponding applicant.

2. Safety note

Warning: This equipment is designed for use in general laboratories, light industrial and office environments. Operation in hazardous areas and/or in a flammable atmosphere is not permitted.

3. General description, application

The Gas Select Module is a PLUGSYS module with a built-in stopcock to select between 4 different gases. The connection to the 4 inlets should be performed with high pressure tubing from a gas tank or wall outlet. The gas pressure must not exceed 8bar (115psi).

With the main lever in the center you can switch between the 4 gasses.

Warning: Never use explosive, flammable or toxic gasses !!!

The outlet of the GSM is connected to the inlet of the PRM (Pressure Regulator Module)




The Pressure Regulator Module PRM is a PLUGSYS module with a built-in pressure regulator. The internal Pressure Regulator is set to regulate the output pressure to 0.5bar (7.2psi)

Warning: Never use explosive, flammable or toxic gasses !!!

The outlet of the PRM must be connected to your setup where you need a low pressurized gas.

4.1 Description GSM Gas Select Module



All input pressure lines should have approximately the same pressure !

If you switch to a non connected input a bit of gas of the pre-selected input may stream out !

4.2 Description PRM Pressure Regulator Module

