# Anesthesia Systems

Complete preconfigured systems

Tabletop or mobile setups

Passive or active scavenging

Customizable to suit your needs







Harvard Apparatus has a proven track record of providing and supporting top quality anesthesia equipment. Anesthesia is a vital part of surgery and we offer a wide array of products to suit your specific needs. We are pleased to offer several complete systems, which are based on the most widely used configurations we see in our customers' laboratories. Our preconfigured packages ensure you have all the parts you need to get started, including vaporizer, regulator, tubing and connectors. If you do not see a combination that is suitable for you, a unique system can be created especially for you.

## Tabletop Systems with Passive Scavenging

- · Ideal flowmeter range for rodents
- · Suitable for mice, rats and other small animals under 10 lb
- Small footprint
- Isoflurane funnel-fill vaporizer with 3-year warranty
- Available with or without induction box
- Passive scavenging



Tabletop Anesthesia System with Induction Box and Passive Scavenging

## Mobile Systems with Passive Scavenging

- · Ideal flowmeter range for rodents
- Suitable for mice, rats and other small animals under 10 lb
- Isoflurane funnel-fill vaporizer with 3-year warranty
- Available with or without induction box
- Passive scavenging
- · Options for multi-animal anesthesia systems



Harvard Apparatus Single Animal Isoflurane System



Harvard Apparatus Mobile Single Animal System



Mobile Anesthesia System with Passive Scavenging (table not included)

## **Tabletop Systems with Active Scavenging**

- · Ideal flowmeter range for rodents
- · Suitable for mice, rats and other small animals under 10 lb
- Small footprint
- Isoflurane funnel-fill vaporizer with 3-year warranty
- Available with or without induction box
- Active scavenging



Tabletop Anesthesia System with MiniVac Active Scavenging

### **Vaporizers**

- Color-coded for isoflurane, halothane and sevoflurane
- Either key-fill or funnel-fill
- · Choice of cage mount or Selectatec mount
- Temperature, flow, and duration of use compensation
- · Capacity of 150 ml



Cage Mount type Vaporizers



Fluosorber Filter Canister



F/Air Filter Canister

### **Passive Scavenging**

- Active charcoal canister packed to optimize the halogenated gas absorbtion
- F/Air can remove up to 50 g and Fluosorber can remove up to 200 g of halogenated gas
- Canister can be thrown away in regular garbage after use



MiniVac Gas Evacuation Unit with F/Air Canister

## MiniVac Active Scavenging

Add the F/Air or Fluosorber to our MiniVac to make an active scavenging system


- · Can be combined with existing setup
- Small footprint
- Adjustable flow
- Quiet
- Easy to use
- Reduces risk of inhalation of residual halogenated gas



MiniVac Gas Evacuation Unit with Fluosorber Canisters



- · Surgical instruments
- Temperature control
- Ventilators
- Lighting
- Magnification
- Operating tables
- Vital signs monitoring
- Infusion products



VentElite Ventilator



Homeothermic Monitoring System



### Contact us for more information!



#### CMA Microdialysis, AB

Torshamnsgatan 30A SE-164 40 Kista, Sweden

phone +46.8.470.10.00 e-mail cma@microdialysis.se web www.microdialysis.se



#### **Harvard Apparatus**

84 October Hill Road Holliston, MA 01746, USA

phone +1.508.893.8999

toll free +1.800.272.2775 (USA Only)

fax +1.508.429.5732 e-mail support@hbiosci.com web www.harvardapparatus.com



#### Harvard Apparatus Canada

6010 Vanden Abeele Saint-Laurent, Quebec H4S 1R9, Canada

phone +1.514.335.0792

toll free +1.800.361.1905 (CAN Only)

fax +1.514.335.3482

e-mail sales@harvardapparatus.ca web www.harvardapparatus.ca



#### **Harvard Apparatus China**

Room 1902E, 19F, Building B Zhong Shan Plaza 1065 West Zhong Shan Road Changning District Shanghai, China

phone +86.21.2230.5128

e-mail china@harvardapparatus.com



#### **Biochrom Limited - Harvard Apparatus UK**

East Wing, Building 1020

Cambourne Business Park, Cambourne Cambridge, CB23 6DW, United Kingdom

phone +44.1732.864001 fax +44.1732.863356

e-mail sales@harvardapparatus.co.uk web www.harvardapparatus.co.uk



#### Hugo Sachs Elektronik / Harvard Apparatus, GmbH

Gruenstrasse 1

March-Hugstetten D-79232, Germany

phone +49.0.7665.9200.0 fax +49.0.7665.9200.90 e-mail info@hugo-sachs.de web www.hugo-sachs.de

## **Panlab**

#### Panlab, S.L. / Harvard Apparatus Spain

C/Energia, 112 08940 Cornellà Barcelona, Spain

phone +34.934.750.697 (International Sales)

+34.934.190.709 (Sales in Spain)

fax +34.934.750.699 e-mail info@panlab.com web www.panlab.com



#### Harvard Apparatus, S.A.R.L.

6 Avenue des Andes Miniparc Building 8 91952 Les Ulis Cedex, France

 phone
 +33.1.64.46.00.85

 fax
 +33.1.64.46.94.38

 e-mail
 info@harvardapparatus.fr

 web
 www.harvardapparatus.fr



